BEMHC MEDICAL PRODUCTS.

Easy-Touch® Glucose Monitoring Systems



Help your customers find their **PERFECT MATCH** with a variety of meters to choose from!



The **EasyTouch**[®] **BluLink[®] meter** is ideal for the user who wants to track and organize their test results. With Bluetooth[®] technology, it can connect right to your smartphone.



The **EasyTouch[®] HealthPro[™] meter** is ideal for the user looking for an upgraded experience with features like a test strip ejector and backlit screen.



The **EasyTouch**[®] **meter** is ideal for those looking for a simple and easy-to-use meter.

Click Here To View APCI Pricing!

CLICK HERE OR SCAN THE QR CODE FOR DIRECT PURCHASING SETUP!





Easy Touch[®] PEN NEEDLES

Comfort You Can Feel!

- Triple-Bevel Cut for a Sharper Point
- Thin Wall Cannulas for an Easier Flow
- Electro-polished to Remove Burrs
- Film-Coated for Maximum Comfort
- Quality Compared to National Leading Brands
- Color-Coded to Gauge
- Easy-to-Read Gauge & Needle Length
- Needles are Made with Surgical Steel

Click Here To View APCI Pricing!

CLICK HERE OR SCAN THE QR CODE FOR DIRECT PURCHASING SETUP!



Easy • Touch[®]

Diabetic Products



MHC Medical Products has been supporting the diabetic community since 2005. Our extensive product offerings and affordable pricing ensure your customers get the supplies they need without breaking the bank! *Call MHC today for details and a complete product list.*

Click Here To View APCI Pricing!

CLICK HERE OR SCAN THE QR CODE



BANC MEDICAL PRODUCTS.

Easy•Touch[®] **FLIPLOCK**[™]

FLIPLOCK NEEDLES 18 1.5" 40mm 50 needles FLIPLOCK NEEDLES Easy • Touch 50 meedles Easy • Touch' 355 FLIPLOCK FL 1″ 25mm 3^{mL}_{CC} SAFFTY SYRINGES Easy Touch 50 Easy Tou F -FLIPLOCK NEEDLES 100 safety syringes 100 200 50 needles

Needle Technology

7 1	
	-

TIGHTEN NEEDLE AND LOAD SYRINGE

Twist the needle to confirm tight fit before use. Pull the plunger back to load the syringe.




INJECT

Push the plunger in one smooth motion to administer medication.



FLIP

Flip the protective shield over the used needle until it clicks into place. Discard according to local regulations



Click Here To View APCI Pricing!

CLICK HERE OR SCAN THE QR CODE FOR DIRECT PURCHASING SETUP!

